

Adapting to Annex 1: Quality Risk Management for Sterile Products

Alan Hoffmeister

Senior Global Scientific Portfolio Specialist

Charles River Laboratories – Microbial Solutions



PDA Aseptic Processing of Biopharmaceuticals Conference 2024

CONNECTING
PEOPLE
AND
SCIENCE
REGULATION®


Contents

1. Annex 1 – Focus on Endotoxin
2. Improved Biofilm Detection to Minimize Impact
3. BET Cartridge Technology Improving Biofilm Detection
4. Cartridge Technology Implementation Customer Case Studies

Annex 1 - Focus on Endotoxin

ANNEX 1 – Critical Utilities

Figure 1 - Filters in a USP Water System



ANNEX 1

MANUFACTURE OF STERILE PRODUCTS

WFI (water for injection) systems are critical utilities for manufacturing pharmaceutical sterile products.

CHAPTER 6 - UTILITIES


Critical utilities are described in point 6.2.:

- Directly contact product.
- Contact materials that will ultimately become part of the product.
- Contact surfaces that come into contact with the product.

RAW MATERIAL

RINSE SOLUTION

WFI as a Critical Utility




Microbiological requirements

Deionized water	
Bioburden limit	NA
Endotoxin limit	NA

Purified water (EP GM0008 / USP 1231)	
Bioburden limit	<100cfu/mL (10000cfu/100ml)
Endotoxin limit	NA

Water for injections (EP GM0169 / USP 1231)	
Bioburden limit	<10cfu/100mL (not sterile)
Endotoxin limit	<0.25EU/mL

Water Quality increase cascade



WFI Microbiological Characteristics

1

Due to lack of nutrients, continuous flow and high temperature (>70°C) in WFI systems not all microorganisms are able to grow in this environment

2

Gram negative rods are more adapted to these environments than other bacteria because they have resistance forms, called BIOFILM

3

As a resistance form, a BIOFILM, has the task to ensure colony survival, is extremely difficult to eliminate and this is why a highly sensitive method for early biofilm detection is critical

4


WFI systems have a bioburden acceptance criteria (less than 10cfu/100ml) this means that these systems are NOT STERILE a certain amount of bioburden is allowed, so BIOFILM formation is a potential (rare, but not impossible) event that has a huge impact to business continuity

Annex 1 – Chapter 6 Utilities – Water systems 6.7.

Due to WFI systems criticality and that WFI requirements accept a certain amount of bioburden coming from the systems that supply water to them, these WFI systems must be:

QUALIFIED

MONITORED



MINIMIZE RISK

microbial contamination/proliferation and pyrogens



PREVENT

formation of **BIOFILMS**

ENSURE

A reliable source of water of an appropriate quality

Microbial Monitoring of CRITICAL UTILITIES - Water Systems

WFI Contamination Control Strategy (CCS) - SAMPLING PLAN

QUALIFICATION

INITIAL

6.8. Water systems should be **qualified** to maintain the appropriate levels microbial control, **taking seasonal variation into account**.

PERIODIC

DAILY

- 6.13. **Regular ongoing microbial monitoring** of water systems should be performed.
- i. All points of use, at a specified interval (QRM).
 - ii. A sample from the point at the end of the distribution loop **each day** that the water is used.

Biofilms in guidelines – EMA Q&A



EUROPEAN MEDICINES AGENCY
SCIENCE MEDICINES HEALTH

1 August 2017
EMA/INS/GMP/443117/2017
GMP/GDP Inspectors Working Group

Questions and answers on production of water for injections by non-distillation methods – reverse osmosis and biofilms and control strategies

Final

PART II - BIOFILMS AND CONTROL STRATEGIES

2. What approach should be taken to maintain control over systems which can be affected by biofilms?

A control strategy should be developed to assess the risks associated with the current manufacturing processes and to determine acceptability of existing control measures. The effectiveness of the **sampling and testing regimes employed at the site should also be critically assessed** in conjunction with the development of a control strategy.

https://www.ema.europa.eu/en/documents/other/questions-answers-production-water-injections-non-distillation-methods-reverse-osmosis-biofilms_en.pdf

Biofilms in guidelines – EMA Q&A

PART II - BIOFILMS AND CONTROL STRATEGIES

3. What is a control strategy in the context of biofilm and contamination control?

A control strategy should take account of the design of the process, the mechanisms required to be put in place to control and ultimately prevent or minimise the risk of contamination.

Such a strategy requires the following thorough process knowledge and understanding taking account of all aspects of contamination control and prevention, including:

- Design
- Water system qualification
- Personnel qualification/training
- Raw Materials, e.g.
- Control strategy including in-process controls applied to
 - Raw Materials - Feed water system - Treatment system
- **Monitoring systems (qualification/calibration) used in the control strategy**

- Preventative maintenance to a standard that will not add significant risk from a contamination view-point
- Robust QMS
 - Deviation handling - Root cause analysis (investigations) – CAPA

Contamination control and steps taken to minimise the risk of contamination are a series of successive linked events/measures. Quality Risk Management tools along with scientific judgement can be applied in determining critical control points.

A contamination control strategy would integrate all of these measures to ensure a more comprehensive approach is taken with respect to prevention and control of microbiological contamination.

Such a strategy should lead to the introduction of a control programme which is an iterative process taking into account all information throughout the lifecycle of the products and processes.

Biofilms in guidelines – EMA Q&A

PART II - BIOFILMS AND CONTROL STRATEGIES

6. Are there any additional measures which should be considered in order to increase the probability of detecting the presence of biofilms?

A robust sampling plan is a requirement. Such a sampling plan forms part of the assessment of the effectiveness of the control strategy employed to minimise such risks of biofilm and general contamination issues. Each potential source of contamination should be incorporated into such a sampling regime. **Ongoing evaluation to determine the appearance of an adverse trend** should be performed, however, the seasonal variation that occurs can only be determined during the annual trend assessment. The effectiveness of an environmental monitoring programme should be formally assessed at minimum on an annual basis.

Sampling programmes for water systems should take account of the quality of the water supply to the system as well as assessing points throughout water generation. Water quality is best assessed through a pre-determined, systematic approach. The loop return should be sampled each day of use of the system in order to provide additional assurance of the quality of water utilised in the manufacturing processes. All points should be sampled on a rotational basis to ensure that the entire system user points are sampled at least once per week.

Routine identification of contaminants isolated during monitoring activities is critical in order to ascertain if there is any shift or change in the flora present within a facility or if certain specific species become more prevalent.

Use of more sensitive endotoxin detection methods should also be taken into account. Alert levels should be set based on the capability of the system and any change or adverse trend should be appropriately investigated.

The frequency of trend analysis and use of trend data is critical. The use of rapid microbiological test methods and systems should be considered in order to improve or increase the probability of early detection and allow timely action to be taken.


Endotoxin testing as a detection method for bacterial biofilms

- American Pharmaceutical Review Article | Dr Tim Sandle

Article Summary: Endotoxin testing can play a role in the earlier detection of biofilms than is possible using conventional bioburden tests. This is on the assumption, albeit one supported by most literature, that much of the bacterial contamination of water systems, and to an extent medical implants, is by Gram-negative bacteria, which on lysis would release endotoxin.

[Endotoxin Testing as a Detection Method for Bacterial Biofilms | American Pharmaceutical Review - The Review of American Pharmaceutical Business & Technology](#)


Improve Biofilm Detection to Minimize Impact



Biofilm impact to production

15 days
to
1 month

COST



GOAL: EARLY BIOFILM DETECTION

To minimize production impact, EARLY BIOFILM DETECTION whilst complying with Annex 1 requirements is needed:

- LESS spread => LESS BIOFILM focus => easier disinfection

- LESS batches impacted



To keep critical utilities under the required quality, CCS requires:

- an effective sampling plan
- and a test method that is specific and sensitive enough to optimize contaminant detection

Methods to assess microbiological quality for WFI systems:



CLASICAL BIOBURDEN METHOD




ENDOTOXIN TESTING AS AN EARLY DETECTION METHOD FOR BACTERIAL BIOFILMS

Methods to assess microbiological quality for WFI systems

	Bioburden	Endotoxin testing
Analyte	Viable microorganisms	Non-viable, highly resistant and stable microorganism particles
Limit of detection	1CFU/200mL	Sensitivity up to 0.005EU/mL parts per trillion (ppt; 10^{-12})
Specificity	All microorganisms that can grow in the media used (R2A?)	High specificity to gram-negative bacteria
Time to result	<ul style="list-style-type: none"> minimum 5 days of incubation ~1-day MANUAL data analysis and report 	<ul style="list-style-type: none"> 15 min to result (CARTRIDGE TECH) MCS (20 samples/hour) Cortex data analysis software

BET Cartridge Technology Improving Biofilm Detection




Improving Biofilm Detection in WFI Systems Using Endosafe® Cartridge Technology

Flexibility

- Random Access
 - Samples can be tested as they arrive to the lab.
 - No need to batch samples > holding time validation

Ease of use

- Easy technician qualification/training
- One error in standard curve preparation for traditional kinetic during a qualification means resampling & retesting of the whole water system



At-line testing

- Remove/reduce sample transportation issues
- Testing next to point of use and data uploaded to Database to be analysed

Sustainability & Compliance

- Method D harmonized method for US and EU
- Uses up to 95% to 100% less horseshoe crab blood derived reagent than traditional kinetic methods

Cartridge technology – How it works

Standard Curve (SC) / Spike Concentration

Negative Product Control (NPC)

Positive Product Control (PPC)

Negative Control (NC)

charles river

Endosafe® - PTS™ Cartridges
Certificate of Analysis

Reorder Code: PTS2001F Cartridge Lot #: 7552158 Expiration Date: Feb 2019

Calibration Code: 114851243599 RSE/CSE Ratio: 13 EU/ug RSE Lot #: H0K354

Archived Standard Curve Range: 1.0,0.1 EU/mL Standard Curve Linearity: 0.998

Standard Curve Mean Reaction Times: 1.0 EU/mL 148 seconds
0.1 EU/mL 402 seconds
0.01 EU/mL 812 seconds

Archived Spike Concentration: 0.110 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES


Qualified Analyst: J. King Date 21 Feb 2017
 Reviewed By: N. Blair Date 27 Feb 2017

Charles River Laboratories, Inc.
1923 Wapoo Road, Suite 43-B
Charleston, SC 29407 USA

CA-PTS20F-04

Duplicate Sample (Well 1 & 3)

Duplicate Sample + endotoxin spike as PPC (Well 2 & 4)



charles river

Endosafe® - PTS™ Cartridges
Certificate of Analysis

Reorder Code: PTS2001F Cartridge Lot #: 7552158 Expiration Date: Feb 2019

Calibration Code: 114851243599 RSE/CSE Ratio: 13 EU/ug RSE Lot #: H0K354

Archived Standard Curve Range: 1.0,0.1 EU/mL Standard Curve Linearity: 0.998

Standard Curve Mean Reaction Times: 1.0 EU/mL 148 seconds
0.1 EU/mL 402 seconds
0.01 EU/mL 812 seconds

Archived Spike Concentration: 0.110 EU/mL Negative Control: Pass

This lot of PTS Cartridges has been tested and meets Quality Control testing requirements for an archived curve, negative controls, and positive product control results.

Store cartridges at 2-25°C. Allow the unopened foil pouch to reach room temperature prior to opening. Cartridges should be used immediately once the foil pouch seal has been opened. Cartridges are for single-test use only.

CAUTION: DO NOT FREEZE THE CARTRIDGES

Qualified Analyst: J. King Date 21 Feb 2017
 Reviewed By: N. Blair Date 27 Feb 2017

Charles River Laboratories, Inc.
1923 Wapoo Road, Suite 43-B
Charleston, SC 29407 USA

CA-PTS20F-04

Advancing **Responsible** Science through Sustainable Endotoxin Testing Solutions

LAL Cartridge Technology



- + **95% Reduction** in horseshoe crab blood from Gel-Clot.
- + FDA-Licensed and compliant with harmonized Pharmacopeia chapters.
- + Robust and scalable technology, from single test to full automation.

rCR Cartridge Technology




- + **100% Replacement** in horseshoe crab blood.
- + Optimized kinetic chromogenic reagent formulated to simulate the natural LAL reaction.
- + Seamless integration with existing suite of scalable, automated, cartridge technology

Regulatory compliance

Complies with USP (<85> <86>) / EP(2.6.14) requirements:

- Kinetic Chromogenic method
- Calibrated to USP RSE, the primary standard
- 2-log ranges with mid-point spikes
- The required sample and spike in duplicates



Customer Case Study 1 – Biotech Efficiency

- During the manufacturing process, both production and QC were involved in the collection of samples for testing
- Often, in-process test results were not available to production within the short time frame required to proceed with downstream processing
- A rapid turnaround of results was required to determine whether endotoxin levels were within limits before proceeding to next step in the manufacturing process
- Six Sigma analysis carried out by objective third party consultant to analyse BET efficiency
- Study analysed sample process flow and turnaround time (TAT) for two sets of in-process buffer samples
- One set tested by central QC using conventional Kinetic Chromogenic plate method
- The other tested using cartridge technology by manufacturing on production floor
- A benefit analysis was performed, taking productivity, quality and financial implications into account

Customer Case Study 1

Method 1




Figure 1a: BET flow chart for in-process buffer samples sent to QC lab using kinetic assays/microplate reader

Method 2




Figure 2a: BET flow chart for in-process buffer samples using the LAL Cartridge/PTS™ at point of use

Customer Case Study 1

Method 1




Figure 1b: TAT for BET using kinetic assay/microplate reader

Method 2




Figure 2b: TAT for BET using LAL cartridges/PTS™

Customer Case Study 1 – Benefit Analysis

Productivity

- Method 2 was shown to reduce total TAT by 1.25 hours per sample set.
- Efficiency was improved by testing at sample collection point.
- Testing performed by manufacturing technicians after a one-day training period.


Quality

- Analyst-to-analyst and lab-to-lab variation was reduced using Method 2 as results were not based on analyst prepared standard curve.
- Retest could be performed immediately eliminating 2.5-hour downtime.
- Simplified OOS investigations due to fewer reagents, steps and accessories required.

Financials

- Employee efficiencies enabled two third-shift employees to transition to other roles.
- One hour reduction in manufacturing cycle time improved capacity and ability to product two additional batches per year.
- Overall efficiency reduced scrap and employee downtime, eliminated waste and decreased cycle time.

Case Study 1: Impact of Cartridge Technology Implementation



Customer Case Study 2

- API produced from Human Urine
- Microbiological “dirty” raw material
 - Exceeds 10^6 cfu/mL
 - High endotoxin levels
- Production process consists of several microbiological and endotoxin reducing steps
- Efficiency of those steps dictates the production process
- Depending on the level of endotoxin contamination rework is needed
- The process continues as production does not have the time to wait for the lab


Case Study 2 - Things to consider....

- *Time to report =*
 - *Time to ship the sample to the lab*
 - *Administrative time*
 - *Sample prep time*
 - *Test incubation time*
 - *Result analysis time*
 - *Verify time*
 - *Approval time*
 - *Time to report the result*
- Product release cycle times are protracted and are the sum of all these sequential activities. That means eight (8) items to work on to speed up the overall testing process

For endotoxin testing this is within one day

Case Study 2 - Process flowchart

- Dedicated to a process & sample type
- Starting material: Human Urine
- Several biochemical column purification
- Samples: Purified protein in different buffers
- Endotoxin concentration must be determined



Case Study 2 - Return on investment?

- RMM PAT Cartridge Technology was implemented
- Cost savings:
 - One batch: \$2 Million
 - Rework: 10% loss in activity (happened in the past)
 - Loss of \$200,000 per batch.
 - 10 batches per year were made.
 - With appropriate in-process analysis rework can be avoided
 - Also testing of the in-coming raw material will avoid bad quality batches
 - No lab needed – all testing performed on production floor

Thank You!

Alan Hoffmeister – Senior Global Scientific Portfolio Specialist
alan.hoffmeister@crl.com