

# Evolution of Aseptic-Containment in Biological and ATMP product filling

**Dr. Hussein Bachir**

**Scientist - GMP Compliance**

**Franz Ziel GmbH**




**2023 Aseptic Processing of Biopharmaceuticals Conference**




# Content

Aseptic-Containment and control measures to achieve OEB 5/6



Strategies for personnel and product safety



Decontamination strategies for aseptic/toxic and biohazard products

Aseptic-Containment and control measures to  
achieve OEB 5/6

## Health Based Exposure Limits

### Toxicological Data

NOAEL

NOEL

### Health Based Exposure Limits

PDE

ADE

### Occupational Limits

OEB

OEL

Health Based Exposure Limits (HBELs) are derived from toxicological data


→ Permitted Daily Exposure (PDE)

$$PDE[\mu\text{g}/\text{day}] = \frac{NO(A)EL\left(\frac{\mu\text{g}}{\text{day}}\right) * \text{Weight Adjustment (kg)}}{F_1 * F_2 * F_3 * F_4 * F_5}$$

Occupational Exposure Limit (OEL): respiratory uptake of permitted daily exposure (PDE)


$$OEL [\mu\text{g}/\text{m}^3] = \frac{PDE}{10 \frac{\text{m}^3}{\text{day}}}$$

## Environmental risks



|                          |                  |   |                             |   |
|--------------------------|------------------|---|-----------------------------|---|
| Very extremely hazardous | OEB 7            | <0.1  | < 10 ng/m <sup>3</sup>      | ACL 7   |
| Extremely hazardous      | OEB 6            | < 1   | < 100 ng/m <sup>3</sup>     | ACL 6   |
| Very highly hazardous    | OEB 5            | 1 – 10  | 100-1000 ng/m <sup>3</sup>  | ACL 5   |
| Highly hazardous         | OEB 4            | 10 – 100  | 1-10 µg/m <sup>3</sup>      | ACL 4   |
| Hazardous                | OEB 3            | 100 – 1000                                      | 10-100 µg/m <sup>3</sup>    | ACL 3   |
| Moderate hazardous       | OEB 2            | 1000 – 10,000                                   | 100-1000 µg/m <sup>3</sup>  | ACL 2   |
| Low hazardous            | OEB 1            | > 10,000  | 1000-5000 µg/m <sup>3</sup> | ACL 1   |
| <b>HAZARD LEVEL</b>      | <b>OEB level</b> | <b>Permitted Daily Exposure (PDEs) (µg/day)</b> |                             | <b>OEL (inhalation) ASEPTIC CONTAINMENT LEVELS (ACLs)</b> |

## Visualisation of OEB exposure amounts



Occupational Exposure Bands (OEBs): Groups of OELs, classification may differ

## primary and secondary containment measures

different systems of containment strategies:

- primary Containments:

physical isolation:


equipment must provide a complete barrier between hazardous process/material and external environment

→ Isolators, bioburden control - API in solution - sterile filtration - positive pressure **vs.** API powder – negative pressure - glove boxes, closed transfer-systems (RTPs), SUS



- secondary Containments:

safeguards and strategies:

support and enhance primary containment measures, particularly in the context of preventing the release of hazardous substances or protecting against contamination (e.g. first air).



## Containment Qualification: Aseptic process filling OEB 5/6 hazardous products





## Strategies for personnel and product safety

Aseptic product  
(toxic)



Personnel Protection  
(Negative Pressure Isolator)



Operator



Aseptic product  
(toxic)



Personnel Protection  
(Negative Pressure Isolator)



GMP

Product Protection  
(Positive Pressure Isolator)


Operator



## Aseptic Containment Strategy (ACS) for operator protection and Cross Contamination Control

- 1 Containment is essential for operator protection and preventing cross-contamination
- 2 Aseptic processing requires an Aseptic-Containment Strategy (ACS) – complementary to CCS
- 3 Containment strategies involve primary and secondary measures working together
- 4 Hazardous product exposure requires characterization with health-based exposure limits (HBELs)
- 5 Aseptic-Containment involves SUS, closed processing, and barrier technology
- 6 Open system isolators use measures like pressure bubbles and directional airflow management for material transfer

## Cross-contamination control



### 1. Product Quality:

crucial to maintain the **purity** and **quality** of products, **contaminants** can be harmful or lead to **product defects**

### 2. Safety:

ensures the **safety** of consumers, workers, and the environment by **preventing** the unintended transfer of **hazardous substances, allergens, or pathogens**

### 3. Regulatory Compliance:

essential for adhering to strict **regulatory standards** and requirements to **protect public health** and ensure **product integrity, guidance of HBELs**

## Viral Containment -Aseptic Containment Strategy (ACS)

**Aseptic-Containment Strategy (ACS) connects to a Contamination Control Strategy (CCS: Annex 1)** when products are biologically active, hazardous or toxic or require containment as a control measure in Cross-contamination control. Both a CCS and ACS need documenting (**critical control points; assessment for monitoring measures**).

**For Viral vector containment a series of primary containment and secondary containment control measures are integrated** to contain the viral vector so it can be decontaminated in-place. Containment measures include: Closed barrier systems, Closed Material transfer systems, Single Use systems and Open Barrier systems with protective airflow 'Bubbles' and airflow patterns.

**In-place decontamination is effected by hydrogen peroxide vapour**, used as a post-production decontamination cycle after batch completion and before open door barrier set up for the next batch.

**Virus may pass through HEPA Filters, so a full 'System' Decon is required** in the post-production decontamination cycle. Barrier system exhaust to outside via double HEPA filtration.

**Viral containment measures apply to airborne contamination (protective airflow patterns) and surface contamination.** Out surfaces of filled units require outer surface Decontamination.


**The Cleanroom, surrounding barrier systems are not a primary Containment** but do provide measures in cross contamination control via PALs/MALs pressure containment 'bubbles'.

# Aseptic/toxic containments are largely, internationally harmonized



Biology ●●●


## BIOSAFETY LEVELs (BSL)



- Foundational principles for biosafety given by WHO Laboratory Biosafety Manual (handling pathogens in research) – level for biosafety and GMP conformity need to be aligned
- Balance risk between operators and patients !
- Biosafety practices often depend on local agencies

## Decontamination strategies for aseptic/toxic and biohazard products






Most Aseptic process filling operations can be considered as localised 'C' contamination for an ACS

# End-point qualification

## Pre-Production-Biodecontamination




- establish Grade A conditions for production process

- End-Point: < 0.1 ppm for biological products

- Short cycle - longer aeration

## Post-Production-Decontamination



- Complete system cycle

- Operator safety & cross-contamination control


- long cycle - short aeration

- End point: < 1 ppm

**Biological products end-point need to be qualified below 0.1 ppm!**

## Exposure studies for end-point validation

- $\text{vH}_2\text{O}_2$  is initially used to achieve Grade A conditions within an Isolator barrier system.
- Pre-Production: end-point qualification, Amplex® Red studies for determination of  $\text{vH}_2\text{O}_2$  residuals
- Post-Production: For viral containment  $\text{vH}_2\text{O}_2$  is used as the principle Post-Production decontamination cycle to facilitate open barrier door access for set-up of next batch.






Figure 1: Mechanism.

- Detection and quantification of H<sub>2</sub>O<sub>2</sub> residues
- HRP catalyzes oxidation of Amplex® Red to highly fluorescent Resorufin
- Emission intensity is proportional to H<sub>2</sub>O<sub>2</sub> concentration

### FZ analytical chemistry laboratory:

- H<sub>2</sub>O<sub>2</sub> measurement device for ppb concentrations (Picarro)
- Amplex red studies for residual analysis can be performed

### On-site:

- after qualification of the gassing cycle, FZ laboratory experts test on site for residuals in operation on the installed Isolator line

### Amplex® Red (10-acetyl-3,7-dihydroxyphenoxazine):

- Fluorogenic probe widely used to detect and quantify H<sub>2</sub>O<sub>2</sub> in biological systems
- MoA: Peroxidase-catalyzed oxidation of Amplex® Red to resorufin.



## Summary

1

Achievement of **OEB 5/6 by Aseptic-Containment Strategy (ACS) that connects to a Contamination Control Strategy (CCS: Annex 1)**

2

For Viral vector containment: combination of **primary and secondary containment** strategies and **qualification is required (PHSS)**

3

Decontamination strategy for **aseptic-toxic API by WIP/CIP**

Decontamination for biological hazards/viral vectors ensured by **Post-Production-Decontamination cycle**

→ Full System Decon is required in Post-Production-Cycle!

4

**Pre-Production:** End-point qualification necessary

Qualified by laboratory analysis!

Thank you !



## Follow up contact details

**Dr. Hussein Bachir - Franz Ziel GmbH  
Scientist - GMP Compliance**

**Hussein.Bachir@ziel-gmbh.com**

