pda.org

The Contribution of Rubber Adhesion to Sealability at Deep Freeze Temperatures

Cathy Zhao, PhD West Pharmaceutical Services, Inc



2020 PDA Asia Pacific Aseptic Processing of Biopharmaceuticals Conference





Purpose

2 Factors impacting vial stopper sealability

Understanding of rubber adhesion at temperatures below glass transition temperature (T_{α})

- Literature search
- Testing
- 4 Takeaways

3

- 5 Future work
- 6 Acknowledgements






Purpose

Elasticity above & below T_g


- What makes rubber seal at temperatures below T_g?
- Shrinkage above & below T_g
- Does polymer chain co-mingling occur between rubber and COC/COP vials?
- The lower the T_g, the better sealability at low temperatures?



COC - cyclic co-polymer; COP - cyclic polymer.



Factors Impacting Vial Stopper Sealability






PDA®



Rubber Stoppers Lose Elasticity Below T_a



DMA – Dynamic Mechanical Analysis



pda.org



Coefficient of Linear Thermal Expansion (CLTE)


<u>SAMPLE</u>	<u>T_g by DMA, °C</u>	<u>CLTE (/°C) x 10⁻⁶,</u> T _g – 125 °C	<u>CLTE (/°C) x 10⁻6,</u> -120 °C - T _g
Rubber 1	-55	293.7	74.1
Rubber 2	-86	305.6	89.87
Rubber 3	-55	221.7	62.84
Rubber 4	-55	222.0	65.10


- Rubber shrinkage varies with formulations.
- The shrinkage is much higher above T_g than below T_g .
- T_g and CLTE do not have direct correlation.






- Physical adsorption intermolecular interactions with intimate contact
 - Wettability
 - COC/COP vs glass with rubber
- Mechanical adhesion interlocking surfaces by filling voids and pores of surfaces
 - Contact area: surface roughness, pressure, viscosity...
- Diffusion merge at interface by diffusion
 - molecules mobile and soluble, oil in thermoplastic elastomer (TPE)
 - polymer chains ends inter-penetration at interface
 - contact time, temperature, molecular weight, physical forms...
- Chemical bonding
- Electrostatic










¹ Adapted from: S Athavale, Adhesion and Must be cited as: Adhesives Theory presentation, September 10, 2010, SIES College of Arts, Science & Commerce, Nerul, Navi Mumbai, India.

pda.org





HNBR - Hydrogenated Nitrile Butadiene Rubber; CS - compression strain; RT - room temperature.




² AG Akulichev, A Tiwari, L Dorogin, AT Echtermeyer and BNJ Persson. Rubber adhesion below the glass transition temperature: Role of frozen-in elastic deformation. EPL, 120 3 (2017) 36002.


pda.org



Role of Adhesion on Rubber Sealability Below $T_{\alpha}^{2,3}$



DMA - Initial compression strain 10% at RT, cooled down to -50°C.



Two different math equations of adhesion force: $T > T_g$ and $T < T_g$


T > T_{g:} JKR (Johnson-
Kendel-Roberts) Theory
$$f_{\rm C} = -\frac{3\pi E_0^* a_{\rm c}^2}{4R}$$
 T < T_{g:} $f_{\rm pull-off} \approx (2\pi\omega a E_1^*)^{1/2}$




² Adapted data from: AG Akulichev, A Tiwari, L Dorogin, AT Echtermeyer and BNJ Persson. Rubber adhesion below the glass transition temperature: Role of frozen-in elastic deformation. EPL, 120 3 (2017) 36002.

³ AG Akulicheva, AT Echtermeyer, BNJ Persson, Interfacial leakage of elastomer seals at low temperatures, International Journal of Pressure Vessels and Piping, 160 2018 (14-23).

Rubber Seal Leak Temperature Determination³ ¹⁰






HNBR: non-filled and 20.4% CB filled Compression: 7% & 30% Sealing surface: milling, grinding & turning $T_{g:}$ ~-16°C (DMA), -23°C (DSC) Vacuum: 100 kPa

Rubber seal leak temperature depends on adhesion force, which depends on:

- rubber formulation
- sealing surface topography and
- compression




³ Data and image adapted from: AG Akulicheva, AT Echtermeyer, BNJ Persson, Interfacial leakage of elastomer seals at low temperatures, International Journal of Pressure Vessels and Piping, 160 2018 (14-23).







Adhesion Force Testing

Adhesion Test Rig







Adhesion Force & Visualization





Rubber Adhesion At Room Temperature With COP & Glass ¹²

Rubber 2 Rubber 1 100 100 Rubber 2, BG, F_{pre} = 5 N Rubber 1, BG, $F_{pre} = 5 \text{ N}$ Rubber 1, BG, $F_{pre} = 12,5 \text{ N}$ 90 90 Rubber 2, BG, F_{pre} = 12,5 N Rubber 1, BG, $F_{pre}^{pre} = 20 \text{ N}$ Rubber 1, BG, $F_{pre}^{} = 120 \text{ N}$ Rubber 2, BG, F_{pre} = 20 N 80 80 Rubber 2, BG, $\epsilon = -30^{\circ}$ Adhesion force Fadh [N] Adhesion force Fadh [N] 70 70 Glass 60 60 50 50 40 Fpre 40 30 30 20 20 10 10 ~ 0 100 10¹ 102 10³ 104 100 10¹ 10² 10³ 104 Contact Time t. [s] Contact Time t, [s] 100 100 Rubber 2, CP, F_{pre} = 5 N ٠ 90 Rubber 1, CP, F_{pre} = 5 N Rubber 2, CP, Fpre = 12,5 N 90 Rubber 1, CP, F = 12,5 N Rubber 2, CP, F_{pre} = 20 N 80 Rubber 1, CP, F_{pre} = 20 N 80 Rubber 2, CP, $\epsilon = -30 \%$ Rubber 1, CP, F Adhesion force Fach [N] 70 COP Adhesion force Fadin [N] 70 60 60 50 50 40 40 30 30 20 20 10 10 0 0 10¹ 10² 100 10¹ 102 10³ 100 10^{3} 104 104 Contact Time t, [s] Contact Time t, [s]


• Adhesion force is rubber formulation-dependent.



PDA®



Deep Freeze Temperature Adhesion Testing Challenges 13








- Equipment limitation
- Time to equilibrium ≥ 4hours

- Ice generation
 - N₂ vapor
 - mold design

- Glue failure
 - steel & COP
 - steel & rubber



Rubber/Glass Adhesion at Low Temperature




- Adhesion force to glass decreases with temperature for both rubbers.
- At 10% 30% compression, no separation between rubber and glass occurred.






Rubber Adhesion to Glass vs COP at Low Temperature ¹⁵







Rubber Adhesion to Glass vs COP at Low Temperature ¹⁶




- Adhesion of Rubber 1 with COP is stronger than with glass at temperature below T_g.
- Rubber 1 generates additional adhesion with COP at temperature below T_g.





Rubber Adhesion to Glass vs COP at Low Temperature



- Adhesion of Rubber 1 with COP is stronger than with glass at temperature below T_g.
- Rubber 1 generates additional adhesion with COP at temperature below T_g.





Takeaways

Adhesion becomes the major sealing force at temperatures below T_g.

Mechanical adhesion generates at T_g, determined by rubber formulation, matrix material, surface roughness and compression pressure.



Rubber shrinkage is much higher at temperatures above $\rm T_g$ than below $\rm T_g.$





Future Work

Test adhesion of CCS system at room temperature vs low temperature.

- Design
- Surface
- Lubrication
- Vacuum

2

Explore math modeling to predict leak temperature.

- Adhesion limit correlation with leak





CCS – Container Closure System





Acknowledgements

Professor Matthias Kröger, Freiberg University of Mining and Technology
Robert Teichert, M.Sc., Freiberg University of Mining and Technology
Le Ho, West Pharmaceutical Services, Inc
Fran DeGrazio, West Pharmaceutical Services, Inc
Allison Radwick, West Pharmaceutical Services, Inc





Thank You!



